

# Biological Water Testing

Reliable safety and quality control

Biological Water Testing



# Ensure Safety and Quality With Reliable and Fast Biological Water Testing

Water is the basis of life. Drinking water, processing and industrial water are relevant resources that need to be controlled and protected carefully.

Comprehensive water analysis is of major importance. Parameters such as total organic carbon (TOC), total nitrogen bound ( $TN_b$ ), and the analysis methods for mercury, trace and toxic elements play an essential role. By supplying high performance TOC/ $TN_b$  analyzers as well as highly sensitive instruments, based on ICP-MS and ICP-OES, Analytik Jena offers a comprehensive solution for chemical water analysis. Protection from infections and unintentional contact with pathogens are also a vital part of water analysis. Reliable and efficient detection methods are essential for optimal action planning and management of hygiene and disinfection.

## Your benefits

- Efficient – from sample to result in a few hours instead of weeks
- Save reaction time – results are available earlier
- Reliable – specific detection kits secure results
- Optimized – perfect combination of device and kit
- Comprehensive – assays for numerous relevant microbiological targets

A conventional cultivating method takes up to 14 days, but only a few hours with Analytik Jena's PCR-based detection method for the highly specific and reliable microbiological parameters (see scheme 1).

### Conservative cultivation



### Molecular biological detection via PCR



## General control parameter: overall bacterial load

The total bacterial load is a key indicator for the hygienic quality of water. It can indicate potential fouling of both, drinking and industrial water.\*

The innuDETECT Bacteria Quantification Assay is a semi-quantitative assay which detects universal bacterial target genes simultaneously with bacterial standard DNA also included in the assay. Internal Controls, co-amplified during the PCR, give maximum assurance on results and validate negative findings.

To demonstrate the excellent quality of the innuDETECT Bacteria Quantification Assay, it was used to analyze the general bacterial load in DNA extracts obtained from a 1 mL water sample. The results are shown in Figure 1. Table 1 demonstrates bacterial loads obtained from real-time PCR, based on exact ct values. Detected bacterial targets in an unknown sample can be quantified upon correlation with standard curves.




Figure 1/Table 1: Co-amplification of bacterial standard DNA (red) allows for semi-quantification of an unknown sample (blue) via ct values. Negative control is shown in black.

| Sample           | Ct value | General bacterial load/reaction |
|------------------|----------|---------------------------------|
| Standard 1       | 21.5     | $1 \times 10^6$                 |
| Standard 2       | 24.8     | $1 \times 10^5$                 |
| Standard 3       | 28.4     | $1 \times 10^4$                 |
| Standard 4       | 32.6     | $1 \times 10^3$                 |
| Standard 5       | 36.1     | $1 \times 10^2$                 |
| Unknown sample   | 31.7     | $2.93 \times 10^3$              |
| Negative control | 40.2     | -                               |

\* Deutsche Trinkwasserverordnung (2001) and 42. Bundesimmissionsschutzverordnung (2017) require determination of culturable microorganisms according to DIN EN ISO 6222:1999-07

### Indicator organisms – a reliable approach for general hygiene control

The presence of specific gastrointestinal microorganisms in water is frequently used as an indicator for fecal contamination. *Escherichia coli*, shiga toxin 1 and/or shiga toxin 2-producing *E. coli* (STEC) as well as shiga toxin-producing *Shigella dysenteriae* can get into water, e.g., via fecal contamination or fertilizers. As this bears the risk of causing severe dysentery, hemorrhagic colitis, and hemolytic uremic syndrome water quality evaluation regarding *E. coli* is obligatory<sup>\*\*</sup>. Several assays, including innuDETECT *E. coli* O157, are designed for the highly specific detection of such pathogens. Figure 2 and Table 2 show a sample application with the combination of three innuDETECT Assays for the characterization of EHEC O104:H4<sup>1</sup>.

<sup>\*\*</sup> Deutsche Trinkwasserverordnung (2001) requires determination of coliforms and *E. coli* according to DIN EN ISO 9308-1:2017-09 & DIN EN ISO 9308-2:2014-06

<sup>1</sup> strain that caused the HUS epidemic in Germany 2011




Figure 3/Table 3: Differential detection of *Legionella* in water samples obtained from cooling towers. *Legionella* spp. (left) are detected in the Cy5 channel and *Legionella pneumophila* (right) in the FAM channel.

|                        | <i>Legionella</i> spp. | <i>Legionella pneumophila</i> |
|------------------------|------------------------|-------------------------------|
| Water sample 1         | (weakly) positive      | negative                      |
| Water sample 2         | negative               | negative                      |
| Water sample 3         | positive               | negative                      |
| Positive control (PCR) | positive               | positive                      |
| Negative control (PCR) | negative               | negative                      |

### Advanced pathogen detection – Ruling out specific hazards

*Clostridia* spp. are bacteria generally found in all environmental habitats including water. Its detection is obligatory according to drinking water regulations<sup>\*\*\*</sup>. The detection via conservative cultivation entails several critical issues. That is why PCR-based assays offer an indispensable alternative. Moreover, *Yersinia enterocolitica* is widely present, especially in animal reservoirs and contaminations are observed in drinking and surface water posing a severe risk to health. Both innuDETECT *Clostridium perfringens* Assay and innuDETECT *Yersinia enterocolitica* Assay offer the perfect solution for fast and sensitive detection of the respective pathogens while the included internal controls ensure result reliability. This is clearly shown in Figure 4 and Table 4.

<sup>\*\*\*</sup> Deutsche Trinkwasserverordnung (2001) requires determination *Clostridium perfringens* according to DIN EN ISO 14189:2016-11




Figure 2/Table 2: Characterization of *E. coli* O104:H4 by real-time PCR via combination of innuDETECT Assays for *E. coli* O104 (red), Shiga Toxin 1 (green), and Shiga Toxin 2 (blue). Shown are signals detected in FAM channel and corresponding results (positive/negative).

| innuDETECT Assay for | Signal FAM channel |
|----------------------|--------------------|
| <i>E. coli</i> O104  | positive           |
| Shiga Toxin 1        | negative           |
| Shiga Toxin 2        | positive           |

### Differential pathogen detection – Co-identification of most relevant subspecies

The presence of *Legionella* needs to be tested in drinking and process water, with a special focus on the subspecies *L. pneumophila*. Legal regulations such as the 42. BImSchV require the evaluation of this Legionnaires' disease-causing pathogen and aims to prevent its emission from evaporation coolers and cooling towers. Identification of *Legionella* by conventional methods is highly challenging and hampered by co-presence of other microorganisms. The innuDETECT *Legionella* Assay overcomes these issues and allows for the simultaneous detection of total *Legionella* spp. and *L. pneumophila* according to the ISO/TS 12857 within a single reaction. Figure 3 and Table 3 show the results of *Legionella* detection and differentiation within water samples from industrial cooling towers. The presence of *L. pneumophila* could only be detected in the positive control of the PCR.




Figure 4/Table 4: Detection signals and ct values in FAM channel during real time PCR were obtained from *C. perfringens* DNA dilutions using the innuDETECT *Clostridium perfringens* Assay in FAM channel.

| DNA amount             | Ct value <i>Clostridium perfringens</i> (FAM) |
|------------------------|---|
| 1 ng                   | 18.44   |
| 100 pg                 | 21.69   |
| 10 pg                  | 25.28   |
| 1 pg                   | 28.99   |
| 100 fg                 | 33.56   |
| Negative control (PCR) | No ct   |

# Order Information

## innuDETECT Assays

| Order number<br>(24 or 96 reactions) | Product name                             | Starting material       | Assay time | qPCR detection channels   | Sensitivity   |
|--------------------------------------|--|-------------------------|------------|---|---|
| 845-IDF-0031024<br>845-IDF-0031096   | innuDETECT Bacteria Quantification Assay | Extracted bacterial DNA | ~ 60 min   | <ul style="list-style-type: none"> <li>▪ FAM (Bacteria)</li> <li>▪ HEX (Internal Control)</li> </ul>  | Dependent on background of negative control   |
| 845-IDF-0033024<br>845-IDF-0033096   | innuDETECT Legionella Assay              | Extracted bacterial DNA | ~ 90 min   | <ul style="list-style-type: none"> <li>▪ FAM (<i>Legionella pneumophila</i>)</li> <li>▪ Cy5 (<i>Legionella spp.</i>)</li> <li>▪ HEX (Internal Control)</li> </ul> | <ul style="list-style-type: none"> <li>▪ <i>Legionella pneumophila</i>: up to 10 DNA copies/PCR</li> <li>▪ <i>Legionella spp.</i>: up to 20 DNA copies/PCR</li> </ul> |
| 845-IDF-0034024<br>845-IDF-0034096   | innuDETECT Clostridium perfringens Assay | Extracted bacterial DNA | ~ 90 min   | <ul style="list-style-type: none"> <li>▪ FAM (<i>Clostridium perfringens</i>)</li> <li>▪ HEX (Internal Control)</li> </ul>  | Up to 10 DNA copies/PCR   |
| 845-IDF-0032024<br>845-IDF-0032096   | innuDETECT Yersinia enterocolitica Assay | Extracted bacterial DNA | ~ 90 min   | <ul style="list-style-type: none"> <li>▪ FAM (<i>Yersinia enterocolitica</i>)</li> <li>▪ HEX (Internal Control)</li> </ul>  | Up to 1 copy/PCR  |
| 845-IDF-0027024<br>845-IDF-0027096   | innuDETECT E. coli O157 Assay            | Extracted bacterial DNA | ~ 60 min   | <ul style="list-style-type: none"> <li>▪ FAM (<i>E. coli</i> O157)</li> <li>▪ HEX (Internal Control)</li> </ul>   | Up to 5 DNA copies/PCR  |
| 845-IDF-0028024<br>845-IDF-0028096   | innuDETECT E. coli O104 Assay            | Extracted bacterial DNA | ~ 60 min   | <ul style="list-style-type: none"> <li>▪ FAM (<i>E. coli</i> O104)</li> <li>▪ HEX (Internal Control)</li> </ul>   | Up to 5 DNA copies/PCR  |
| 845-IDF-0025024<br>845-IDF-0025096   | innuDETECT Shiga Toxin 1 Assay           | Extracted bacterial DNA | ~ 60 min   | <ul style="list-style-type: none"> <li>▪ FAM (Shiga toxin 1)</li> <li>▪ HEX (Internal Control)</li> </ul>   | Up to 5 DNA copies/PCR  |
| 845-IDF-0026024<br>845-IDF-0026096   | innuDETECT Shiga Toxin 2 Assay           | Extracted bacterial DNA | ~ 60 min   | <ul style="list-style-type: none"> <li>▪ FAM (Shiga toxin 2)</li> <li>▪ HEX (Internal Control)</li> </ul>   | Up to 5 DNA copies/PCR  |

### Headquarters

Analytik Jena AG  
Konrad-Zuse-Str. 1  
07745 Jena · Germany

Phone +49 3641 77 70  
Fax +49 3641 77 9279  
info@analytik-jena.com  
www.analytik-jena.com

Pictures: Analytik Jena AG, © pixabay - Pexels  
Subjects to changes in design and scope of delivery as well as further technical development!