# Bedienungsanleitung/Operating Manual Autoinjector (Typ/type AI-EA)

---

Hersteller                   Analytik Jena GmbH+Co. KG  
Konrad-Zuse-Str.1  
07745 Jena · Deutschland  
Telefon + 49 3641 77 70  
Fax       + 49 3641 77 9279  
E-Mail   info@analytik-jena.com

Service                     Analytik Jena GmbH+Co. KG  
Konrad-Zuse-Str. 1  
07745 Jena · Deutschland  
Telefon + 49 3641 77 7407  
Fax       + 49 3641 77 7449  
E-Mail   service@analytik-jena.com

Allgemeine                <http://www.analytik-jena.com>  
Informationen

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General Information   <http://www.analytik-jena.com>

Edition                B (07/2023)

Technical documentation  
made by               Analytik Jena GmbH+Co. KG

© Copyright 2023, Analytik Jena GmbH+Co. KG

# Contents

|           |                                                                  |           |
|-----------|------------------------------------------------------------------|-----------|
| <b>1</b>  | <b>Basic information .....</b>                                   | <b>5</b>  |
| 1.1       | Notes on these operating instructions .....                      | 5         |
| 1.2       | Intended use .....                                               | 6         |
| <b>2</b>  | <b>Safety instructions.....</b>                                  | <b>7</b>  |
| 2.1       | General notes.....                                               | 7         |
| 2.2       | Safety markings on the autoinjector .....                        | 7         |
| 2.3       | Requirements for the operating personnel.....                    | 7         |
| 2.4       | Safety instructions, transport and commissioning .....           | 7         |
| 2.5       | Safety instructions for operation .....                          | 8         |
| 2.5.1     | General .....                                                    | 8         |
| 2.5.2     | Safety instructions – protection against explosion and fire..... | 8         |
| 2.5.3     | Safety instructions – electrical equipment.....                  | 8         |
| 2.5.4     | Handling of auxiliary and operating materials and samples .....  | 8         |
| 2.5.5     | Safety instructions – maintenance and repair .....               | 8         |
| 2.6       | Safety equipment / behavior during emergencies .....             | 9         |
| <b>3</b>  | <b>Function and setup.....</b>                                   | <b>10</b> |
| <b>4</b>  | <b>Installation and commissioning .....</b>                      | <b>13</b> |
| 4.1       | Installation conditions.....                                     | 13        |
| 4.1.1     | Environmental conditions.....                                    | 13        |
| 4.1.2     | Energy supply.....                                               | 13        |
| 4.1.3     | Device layout and space requirements.....                        | 13        |
| 4.2       | Connecting the autoinjector to the elemental analyzer.....       | 13        |
| <b>5</b>  | <b>Operation.....</b>                                            | <b>15</b> |
| <b>6</b>  | <b>Maintenance and care.....</b>                                 | <b>18</b> |
| 6.1       | Adjusting the autoinjector for the syringe .....                 | 18        |
| <b>7</b>  | <b>Troubleshooting .....</b>                                     | <b>19</b> |
| 7.1       | Displays in the Device Status window .....                       | 19        |
| 7.2       | Equipment faults and analytical problems.....                    | 19        |
| <b>8</b>  | <b>Transport and storage.....</b>                                | <b>20</b> |
| 8.1       | Preparing the device for transport .....                         | 20        |
| 8.1.1     | Packing the autoinjector .....                                   | 20        |
| 8.2       | Storage .....                                                    | 21        |
| <b>9</b>  | <b>Disposal .....</b>                                            | <b>22</b> |
| <b>10</b> | <b>Specification .....</b>                                       | <b>23</b> |
| 10.1      | Technical data.....                                              | 23        |
| 10.2      | Guidelines and standards.....                                    | 24        |

# Figures

Fig. 1 multi EA 5100 with autoinjector ..... 10

Fig. 2 Autoinjector setup ..... 11

Fig. 3 Base plate (from below) ..... 11

Fig. 4 Microliter syringe (100 µl)..... 12

Fig. 5 Autoinjector connection to the elemental analyzer ..... 12

Fig. 6 AI-EA Adjustment window..... 18

# 1 Basic information

## 1.1 Notes on these operating instructions

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Contents                 | <p>The AI-EA autoinjector is a system module of the modular elemental analyzers multi EA 5000 and multi EA 5100. This user manual is therefore only valid in conjunction with the following documents:</p> <ul style="list-style-type: none"> <li>▪ User manual for the multi EA 5000 and multi EA 5100</li> <li>▪ Software manual for the multiWin control and analysis software</li> </ul> <p>These operating instructions provide information about the design and operation of the sample feeding module and also provide operating personnel with the necessary know-how for the safe handling of the device and its components. The operating instructions also provide information on maintenance and service, as well as on possible causes and remedies in the event of malfunctions.</p> |
| User requirements        | <p>These operating instructions are intended for users who are familiar with the principles of elemental analysis. Knowledge of the safety guidelines for working in a chemical laboratory and for handling the chemicals that are used is required. Furthermore, the user must have basic knowledge of working with a computer.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Conventions              | <p><b>Instructions for actions</b> which occur in chronological order are numbered and combined in action units.</p> <p><b>Warnings</b> are indicated by warning triangles and a signal word. The type, source and consequences of the danger are stated together with notes on preventing the danger.</p> <p>The elements of the control and analysis program are indicated as follows:</p> <ul style="list-style-type: none"> <li>▪ Program terms are indicated by small caps (e.g., SYSTEM menu).</li> <li>▪ Buttons are shown by square brackets (e.g., [OK]).</li> <li>▪ Menu items are separated by arrows (e.g., SYSTEM ► DEVICE).</li> </ul>                                                                                                                                               |
| Symbols and signal words | <p>The operating instructions use the following symbols and signal words to indicate hazards or instructions. The warnings are always placed before an action.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

### WARNING

Indicates a potentially hazardous situation which might cause death or very serious injuries (deformities).

### CAUTION

Indicates a potentially hazardous situation which might cause light or minor injuries.

### Attention

Provides information on potential material or environmental damage.

## 1.2 Intended use

The autoinjector is a sample feeding module of the modular elemental analyzers and multi EA 5000 and multi EA 5100. The autoinjector injects individual liquid samples directly into the combustion tube. It can be used only with the analyzer in vertical operation.

The autoinjector may be used only for the methods for the determination of sulfur, nitrogen, chlorine and carbon content in liquid samples described in these operating instructions and in the operating instructions for the multi EA 5000 and multi EA 5100. Any other use is not as intended! The operator is solely responsible and liable for any damages that result from such unintended use.

The following substances must not be dosed into the combustion tube of the multi EA 5000 and multi EA 5100 with the autoinjector because there is a risk of explosion:

- highly flammable organic compounds (e.g., isopentane)
- substances tending to spontaneous decomposition (e.g., peroxides)
- explosives, explosive materials (e.g., trinitrotoluene, inorganic azides)

## 2 Safety instructions

### 2.1 General notes

For your own safety and to ensure error-free operation of the device, please read this chapter carefully before commissioning.

Observe all the safety instructions provided in these operating instructions and in the operating instructions for the multi EA 5000 and multi EA 5100 as well as all messages and instructions displayed by the multiWin control and analysis program on the monitor.

### 2.2 Safety markings on the autoinjector

Notice symbols have been attached to the autoinjector. These must be strictly observed.

Damaged or missing symbols can cause incorrect actions leading to personal injury or material damage! Symbol labels must not be removed or wetted with organic solvents! Damaged symbol labels must be replaced without delay!

| Symbol                                                                              | Meaning                                                                                                                                                                                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Risk of burns at the hot furnace!<br>Make sure not to touch the hot combustion system when installing or removing the autoinjector.                                                                                                                                 |
|  | The device contains controlled substances. Analytik Jena guarantees that, if the device is used as intended, these substances will not escape from the device in the next 25 years and will not pose a threat to the environment or health within this time period. |

### 2.3 Requirements for the operating personnel

The autoinjector may only be operated by qualified specialist personnel instructed in the use of the device. The instruction also includes conveying the contents of these operating instructions. We recommend training by qualified employees of Analytik Jena or its representatives. The operating instructions must be accessible to the operating and service personnel.

### 2.4 Safety instructions, transport and commissioning

Observe the following notes:

- Insufficiently secured components pose a risk of injury! During transport the components of the equipment must be secured in accordance with the specifications in these operating instructions.

## 2.5 Safety instructions for operation

### 2.5.1 General

The operator of the autoinjector must make sure before each commissioning that the condition of the device is sound. This applies in particular after each modification or extension of the device or its repair.

Observe the following notes:

- The device may only be operated if all items of protective equipment (covers in front of electronic components) are in place, properly installed and fully operational.
- The needle of the microliter syringe has a very sharp tip. Do not touch the syringe by the needle. This could result in injuries or contamination of the needle.
- The furnace of the elemental analyzer operates with temperatures of up to 1100 °C. The hot syringe needle must not be touched during or immediately after the injection.
- Ensure that no liquid reaches the plug connectors or enters the interior of the device! Danger of electric shock!
- Caution when handling glass components. There is the risk that glass components may break and cause injuries!
- Keep all combustible materials away from the device.

### 2.5.2 Safety instructions – protection against explosion and fire

The autoinjector must not be operated in an explosive environment. Smoking or handling open flames are prohibited in the room in which the device is operated!

### 2.5.3 Safety instructions – electrical equipment

Observe the following notes:

- The connection cable between the autoinjector and the elemental analyzer may only be connected or disconnected when the device is switched off. Connect the autoinjector only via the provided interface to the elemental analyzer.

### 2.5.4 Handling of auxiliary and operating materials and samples

The operator is responsible for the selection of substances used in the process as well as for their safe handling. This is particularly important for radioactive, infectious, poisonous, corrosive, combustible, explosive and otherwise dangerous substances.

### 2.5.5 Safety instructions – maintenance and repair

Observe the following notes:

- The operator may only carry out the tasks listed in chapter "Maintenance and care" p. 18.

- The exterior of the autoinjector may only be cleaned using a cloth slightly moistened with water, but not dripping wet.
- Any maintenance on the device may only be carried out when the device is disconnected from the power supply (unless stated otherwise).
- Only accessories and original replacement parts from Analytik Jena may be used.

## 2.6 Safety equipment / behavior during emergencies

Observe the following notes:

- During dangerous situations or after accidents have occurred, if there is no immediate risk of injury, where possible immediately switch off the elemental analyzer at the power switch and/or disconnect the power plug from the power socket!
- Close the gas supply as immediately as possible after switching off the device!

### 3 Function and setup

The autoinjector injects individual samples with a volume of 1 to 100 µl directly into the combustion tube of the elemental analyzer.

The samples are fed semi-automatically:


- The user draws up the liquid sample into the supplied microliter syringe.
- The user then inserts the filled syringe into the autoinjector.
- The autoinjector doses the sample into the combustion tube very precisely and reproducibly at the preset dosing speed.

The autoinjector is inserted into the opening at the top of the elemental analyzer. It can be used only with the analyzer in vertical operation. The autoinjector can be optionally fitted with microliter syringes with a volume of 50 µl or 100 µl. The sample feeding module is controlled via the multiWin software.

#### Setup


Fig. 1 multi EA 5100 with autoinjector

**Fig. 2** Autoinjector setup

- |   |                                                                              |   |                  |
|---|------------------------------------------------------------------------------|---|------------------|
| 1 | Base element with mechanics and control electronics                          | 5 | Syringe bracket  |
| 2 | Guiding device to push down the piston of the syringe in a controlled manner | 6 | Syringe cylinder |
| 3 | Syringe piston                                                               | 7 | Needle guide     |
| 4 | Conical guide at the syringe with positioning pin                            | 8 | Base plate       |

The autoinjector can be moved on the base plate until its guide is positioned exactly above the injection port (septum with screw cap) of the elemental analyzer.

**Fig. 3** Base plate (from below)

- |   |                                                     |
|---|-----------------------------------------------------|
| 1 | Guide for the injection port                        |
| 2 | Adjustable mounting for the AI-EA on the base plate |

The autoinjector may only be fitted with the microliter syringes provided for this purpose by Analytik Jena. The syringes have a conical guide with a positioning pin that fits into the syringe bracket of the autoinjector. The needle can be replaced if, for example, it is blocked.


Fig. 4 Microliter syringe (100 µl)

- |   |                          |   |                                    |
|---|--------------------------|---|------------------------------------|
| 1 | Exchangeable needle      | 3 | Conical guide with positioning pin |
| 2 | Graduated glass cylinder | 4 | Piston                             |

The elemental analyzer controls the autoinjector via the interface on the rear of the device. At the same time the elemental analyzer supplies the autoinjector with an operating voltage of 24 V.

Fig. 5 Autoinjector connection to the elemental analyzer

## 4 Installation and commissioning

### 4.1 Installation conditions

#### 4.1.1 Environmental conditions

The location requirements correspond to those of the multi EA 5000 and multi EA 5100 analyzer (→ multi EA 5000/EA 5100 operating instructions).

#### 4.1.2 Energy supply

The elemental analyzer supplies the autoinjector with a voltage of 24 V. Connect the autoinjector only via the provided interface to the elemental analyzer.

#### 4.1.3 Device layout and space requirements

The autoinjector is mounted at the top of the elemental analyzer. The required height results from the height of the analyzer and the height of the sample feeding module. There must be a clearance of at least 10 cm between the device system and any shelf or cabinet above the system.

| Device                               | Dimensions (width x height x depth) in cm |
|--------------------------------------|-------------------------------------------|
| multi EA 5000/EA 5100 (basic device) | 53 x 47 x 56                              |
| Autoinjector                         | 15 x 27 x 24                              |

### 4.2 Connecting the autoinjector to the elemental analyzer

#### CAUTION

Risk of burns at the hot furnace! Switch off the analyzer and let the analyzer cool down before installing the autoinjector.

#### Attention

The autoinjector may only be connected after the analyzer has been switched off at the power switch and disconnected from the power supply!

Do not hold the autoinjector by the red front panel when attaching it or removing it.

1. Prepare the analyzer for vertical operation.
2. Screw the base plate to the autoinjector using the 2 hexagon socket screws (→ 2 in Fig. 3 p. 11).
3. Insert the autoinjector into the opening at the top of the analyzer. The autoinjector is correctly inserted when the base plate completely sits in the opening and is flush fitting with the surface of the analyzer.

- Establish the electrical connection between the autoinjector and the analyzer.

- ✓ This completes the installation of the autoinjector.

## 5 Operation

### CAUTION

Risk of injury when handling the needle

The needle of the microliter syringe has a very sharp tip. Do not touch the syringe by the needle. This could result in injuries or contamination of the needle.

### Attention

The autoinjector may only be fitted with the syringes approved for this purpose.

When inserting the syringe into the autoinjector, make sure that you do not touch the piston. Otherwise dosing could start prematurely.

For a detailed description of the multiWin control and evaluation software, see the multiWin operating instructions.

1. Switch on the PC. Establish the gas supply for the elemental analyzer.
2. Switch on the elemental analyzer and the detector at the power switch.
  - ✓ The LED on the left front door of the analyzer illuminates green after approx. 30 s.
3. Start the multiWin control and evaluation software on the PC.
4. Create a method via the METHOD ► METHOD – NEW menu item or activate a method via the METHOD ► METHOD – ACTIVATE menu item.
5. Specify the sample volume in the method on the DETERMINATIONS tab under OPERANDS ► SAMPLE VOLUME (range 1 to 100 µl).
6. Set the dosing speed on the tab PROCESS/SYSTEM tab under SAMPLES: DOSE (range 0.2 to 2 µl/s).

**i** Attention! The following method settings are not relevant for the autoinjector. They only apply to the MMS 5000/MMS 5100 autosampler.

- Settings on the DETERMINATIONS tab under RINSE
- The SAMPLE: DRAW UP setting on the PROCESS/SYSTEM tab
- Settings on the PROCESS/SYRINGE tab

7. Start the analysis by clicking on [START MEASUREMENT] on the main window.
8. Create a new analysis group or select an existing analysis group. Confirm with [OK].
9. Create an analysis sequence.

When using the autoinjector, the MANUAL MEASUREMENT option is preset in the sequence under RACK/SELECTION.

10. Enter the sample ID in the NAME field.
11. Complete the lines of the sequence in accordance with the sampling process.

If the autoinjector should dose a different sample volume than specified in the method: Enter the desired dosing volume between 1 µl and 100 µl under SAMPLE

VOLUME. In the sequence, a different dosing volume can be specified for each sample.

12. Click on  in the sequence row to release a row for analysis. Alternatively: Release all sequence rows by clicking on .

After releasing the sequence entries, click on [OK].


13. Start the measurement: In the MEASUREMENT – ANALYSIS window, click on [START MEASUREMENT].

The software prompts you to insert the syringe.

14. Draw the sample into the syringe of the autoinjector without bubbles.

For this purpose, immerse the needle into the sample solution, draw up and eject sample solution several times with high speed in order to get rid of tiny bubbles. Then draw up the required volume slowly. Carefully wipe the outside of the needle with cellulose. Otherwise, sample residues could contaminate the septum.

15. Remove the red cover from the autoinjector.

**i** Attention! When inserting the syringe into the autoinjector, do not touch the piston and do not press it down!

16. Hold the syringe by the white plastic cone (2) and guide it through the needle guide (4) of the autoinjector. Push the needle through the septum of the combustion tube.

17. Insert the needle into the combustion tube until the cone of the syringe fully rests on the holder of the autoinjector (3).

The black positioning pin must be located in the recess of the holder.

18. Reattach the red cover.

19. In the dialog window, confirm the insertion of the syringe with [OK].

- ✓ The autoinjector moves the piston (1) down under software control and injects the samples into the combustion tube. Analysis is carried out.

20. When prompted by the software, immediately remove the syringe from the autoinjector. Prepare the next measurement.

---

### Attention

Risk of erroneous results! Make sure that you draw up the sample volume specified in the method or sequence in the syringe.

The autoinjector is equipped with a sensor that detects the position of the piston. It injects the actually drawn sample volume into the combustion tube. However, the multiWin software calculates the element content of the sample based on the sample volume specified in the method/sequence.

When prompted by the software, immediately insert the syringe into the autoinjector or remove it again. Irregularities in handling will increase the scattering of the measurement results. In addition, the septum and the needle will age quickly if the syringe remains inserted for too long.

---

## 6 Maintenance and care

| Maintenance task                      | Maintenance interval                                   |
|---------------------------------------|--------------------------------------------------------|
| Clean and maintain the device         | Weekly                                                 |
| Clean the needle with a cleaning wire | For example, if clogged with particles from the septum |
| Replace the needle                    | If damaged or if the blockage cannot be removed        |
| Replace the syringe                   | As required                                            |

### 6.1 Adjusting the autoinjector for the syringe

During the software-supported adjustment of the autoinjector, the system checks the geometry of the syringe that is used. If required, the number of steps that the stepper motor moves down the syringe piston is adjusted to the geometry of the syringe.

100 µl microliter syringe

1. Select the SYSTEM ► AI-EA ADJUSTMENT command.


Fig. 6 AI-EA Adjustment window

- In the AI-EA ADJUSTMENT window, select [ADJUSTMENT CHECK] to check the adjustment.
  - Select [SYRINGE ADJUST] to start the adjustment of the syringe.
2. When prompted, pull out the syringe piston exactly to the maximum volume of 100 µl. Insert the syringe into the autoinjector and confirm with [OK].
    - ✓ The autoinjector pushes the piston down to half the sample volume.
  3. Read the actual value (e.g., 48 µl) from the graduated syringe and enter it in the input field. Confirm with [OK].
    - ✓ The autoinjector has been adjusted.
  4. Select [ADJUSTMENT CHECK] to test whether the 50 µl mark is exactly reached during the next dosing.

50 µl microliter syringe

Draw up the syringe plunger to 50 µl during adjustment. Multiply the actual value read from the syringe cylinder by a factor of 2 and enter it in the input field.

Repeat the adjusting with each syringe change.

## 7 Troubleshooting

### 7.1 Displays in the Device Status window

During initialization, the multiWin control software recognizes the connected sample feeding modules and displays them in the STATUS ANALYZER window. The autoinjector is recorded with its syringe size: Display: AI-EA SYRINGE: 100 µL.

The following status display indicates malfunctions:

| Display                                            | Description                                                                                                                                                       |
|----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No sample feeding module displayed (display empty) | Device not ready for measurement, autoinjector not recognized: <ul style="list-style-type: none"> <li>■ Connect the autoinjector</li> <li>■ Initialize</li> </ul> |

### 7.2 Equipment faults and analytical problems

Other problems that are not detected by system monitoring can also occur. Starting a measurement is possible. Such errors are usually recognized by analytical problems.

If the suggested solutions are not successful, inform the customer service department of Analytik Jena.

| Error                      | Possible cause                                                  | Remedy                                                                   |
|----------------------------|-----------------------------------------------------------------|--------------------------------------------------------------------------|
| Lower results              | Dosing fault                                                    | Check the filling of syringe, adjustment and dosing                      |
| Lower / higher results     | Sample volume deviates from the volume settings in the software | Adjust sample volume to the software settings                            |
| Carryover                  | Inadequate rinsing of the dosing syringe                        | Rinse the dosing syringes adequately prior to sampling                   |
|                            | Dosing faulty                                                   | Check dosing                                                             |
|                            | Contamination in the combustion system                          | Check the combustion system and clean it if necessary                    |
| Scattering measured values | Dosing faulty                                                   | Check dosing                                                             |
|                            | Clogged needle                                                  | Use the wire provided for this purpose to clean the needle or replace it |

## 8 Transport and storage

### 8.1 Preparing the device for transport

---

#### WARNING

Risk of damage to health due to improper decontamination!

Perform a professional and documented decontamination of the device before returning it to Analytik Jena. The decontamination report is available from the customer service department when registering the return. Analytik Jena must refuse acceptance of contaminated devices. The sender may be liable for any damage caused by inadequate decontamination of the device.

---

#### CAUTION

Risk of burns at the hot furnace! Make sure not to touch the hot combustion system when removing the autoinjector.

When removing the glass components there is a risk of injury from glass breakage! Carefully remove the microliter syringe from the analyzer!

Risk of injury when handling the needle! The needle of the syringe has a very sharp tip. Do not touch the syringe by the needle when removing it.

---

#### 8.1.1 Packing the autoinjector

---

#### Attention

The autoinjector may only be disassembled after the analyzer has been switched off at the power switch and the power plug has been disconnected from the power socket!

Do not hold the autoinjector by the red front panel when attaching it or removing it.

---

1. Disconnect the electrical connection on the rear of the analyzer.
2. Remove the microliter syringe from the autoinjector and clean it. Pack the syringe separately in the packaging supplied.
3. Remove the autoinjector from the elemental analyzer.
4. Detach the base plate via the 2 hexagon socket screws on the underside.
5. Pack the autoinjector and the base plate in their original packaging.

## 8.2 Storage

---

### Attention

Environmental influences and condensate formation can destroy individual components of the device!

The device must only be stored in air-conditioned rooms. The atmosphere must be low in dust and free from aggressive vapors.

---

If the device is not installed immediately after delivery or not required for prolonged periods, it should be stored in its original packaging. A suitable desiccant should be added to the equipment to prevent damage from moisture.

Ambient conditions

For information on the requirements with regard to the ambient conditions at the storage location, refer to "Technical data" p. 23.

## 9 Disposal

At the end of its service life, the autoinjector and all its electronic components must be disposed of as electronic waste in accordance with the applicable regulations.

## 10 Specification

### 10.1 Technical data

|                                     |                                     |                            |
|-------------------------------------|-------------------------------------|----------------------------|
| Ambient conditions during operation | Temperature range                   | 21 to 35 °C                |
|                                     | Humidity                            | Max. 90% at +30 °C         |
|                                     | Air pressure                        | 0.7 to 1.06 bar            |
| Ambient conditions during storage   | Temperature range                   | 15 to 55 °C                |
|                                     | Humidity                            | 10 to 30 % (use dessicant) |
| General characteristics             | Designation/type                    | Autoinjector (AI-EA)       |
|                                     | Dimensions (width x height x depth) | 15 x 27 x 24 cm            |
|                                     | Mass                                | 1.5 kg                     |
|                                     | Syringe                             | 100 µl or optionally 50 µl |
|                                     | Sample volume                       | 1 to 100 µl                |
| Electrical variables                | Connection                          | with 24 V = 2.0 A          |

## 10.2 Guidelines and standards

|                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Protection class, protection type | The autoinjector belongs in protection class III. The housing has safety type IP 20.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Device safety                     | <p>The autoinjector complies with the safety standards</p> <ul style="list-style-type: none"><li>▪ EN 61010-1</li><li>▪ EN 61010-2-081</li></ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| EMC compatibility                 | <p>The autoinjector has been checked for interference emission and noise immunity.</p> <ul style="list-style-type: none"><li>▪ It meets the requirement for interference emissions according to EN 61326-1 (EN 55011 group 1, class B)</li><li>▪ It meets the requirements for noise immunity according to EN 61326-1 (requirements for use in a basic EMC environment)</li></ul>                                                                                                                                                                                                                                 |
| EU directives                     | <p>The autoinjector is built and tested according to standards which comply with the requirements of the EU Directives 2014/30/EU (EMV-RL) and 2011/65/EU (RoHS). The device leaves the factory in a sound condition as far as technical safety is concerned. To maintain this condition and to ensure safe operation, the operator must strictly observe the safety and operating instructions contained in this operating manual. For accessories delivered with the device and system components from other manufacturers, the information provided in the respective operating instructions has priority.</p> |
| Guidelines for China              | <p>The device contains substances subject to regulation (according to directive "Management Methods for the Restriction of the Use of Hazardous Substances in Electrical and Electronic Products"). Analytik Jena guarantees that, if the device is used as intended, these substances will not leak within the next 25 years and will not pose a threat to the environment or health within this time period.</p>                                                                                                                                                                                                |