21 CFR Part 11 Compliance mit qPCRsoft

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


Role	Permissions
☒ Administrator	☒ User is allowed to manage settings
☒ Supervisor	☒ User is allowed to manage export configurations
☒ Operator	☒ User can configure devices
☒ Data auditor	☒ User is able to edit, create and delete other users
☒ Service engineer	☒ User can edit color modules
	☐ User can manage data groups
	☒ User can change the licence of the software
	☐ User can import LIMS
	☐ User can stop and start measurement
	☐ User can calculate Ct/Tm
	☐ User can use unrelease methods for measurement
	☐ User can manage templates
	☐ User can edit template general parameters
	☐ User can edit template thermal cycler parameters
	☐ User can edit template scan parameters
	☐ User can manage projects
	☐ User can add and edit analyses
	☐ User can import and export analyses

Benutzermanagement

Created	User	Text	Group	Category
08.08.2025 10:15:22	Administrator	Administrator recovery generated...	User Audit Trail	Login
08.08.2025 10:15:22	Administrator	Administrator recovery generated...	User Audit Trail	Login
08.08.2025 10:15:15	Phantomas (Phantom)	User Admin Administrator recovery generate...	User Audit Trail	Login
08.08.2025 10:14:12	Phantomas (Phantom)	Color module config. changed for Virtual iris...	Device Audit Trail	Operate
08.08.2025 10:12:46	Phantomas (Phantom)	Color module config. changed for Virtual iris...	Device Audit Trail	Operate
08.08.2025 10:09:02	Phantomas (Phantom)	Project PowerQuant opened	Project Audit Trail	Result
08.08.2025 10:08:48	Phantomas (Phantom)	Device parameters changed for during the d...	Device Audit Trail	Service
08.08.2025 10:01:58	Phantomas (Phantom)	Project PowerQuant opened	Project Audit Trail	Result
08.08.2025 10:01:40	Phantomas (Phantom)	Project Absolute Quantification created	Project Audit Trail	Result
08.08.2025 09:33:47	Phantomas (Phantom)	Device parameters changed for during the d...	Device Audit Trail	Service
08.08.2025 09:33:36	Phantomas (Phantom)	Device Virtual iris device of type 96 connecte...	Device Audit Trail	Operate
08.08.2025 09:24:32	Phantomas (Phantom)	Instrument 3107B-0759 initialized	Device Audit Trail	Operate
08.08.2025 09:24:32	Phantomas (Phantom)	Device 3107B-0759 connected via usb://\./u...	Device Audit Trail	Operate

Audit-Trail

Please sign the action:

Please confirm your credentials to sign:

Username: Testuser Password: [input field]

Version: [input field] Is Current Version

Sign Status: Signatures

Template: [Created] [Checked] [Released] [Invalid]

Project: [Created] [Checked] [Released] [Invalid]

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© Analytik Jena GmbH+Co. KG
Analytik Jena GmbH+Co. KG
Bilder: Analytik Jena GmbH+Co. KG