

We Normalize Speed

CyBio FeliX SELECT Head





The SELECT Head revolutionizes modular, benchtop liquid handling with its 8 independent channels and volumes. Simplify complex transfers by executing directly from your CSV files, liberating you from complicated programming. Avoid slow pipetting bottlenecks in your workflows and automate faster than with a single tip.


What you get:

- **Increased efficiency:** Tasks can be completed faster and with less effort.
- **Pre-prepared applications:** Workflows include basic data handling and reporting, but options such as custom handling and calculations are possible.
- **Increased adaptability:** The variable channels accommodate different liquid volumes, simplifying sample handling, while the SELECT Head can quickly rearray into various patterns or plate types.

Key Supported Applications


Hit Picking just got faster

Select any 8 wells simultaneously across 96 and 384 formats. Speed up your workflow by reducing repetitive movements.



Normalization just got faster

Normalize samples quickly to the correct concentration and volume. The SELECT Head can aspirate up to 8 different volumes of sample and diluent in a single transfer. Automatically determine required volumes from your input files, eliminating manual calculations.



Need help getting started?

Our motivated application team is here to support you with any challenges you face. The flexible SELECT Head can handle a wide range of tasks, ensuring we always find a solution for your needs. Don't hesitate to reach out and let us help you achieve your goals.

Contact sales@analytik-jena.com

What is Your Challenge?

Choose the right SELECT Head

With volume ranges from 1 µL to 1000 µL, there is a SELECT Head to automate every application. Combine the SELECT Head with any new or existing CyBio FeliX instrument, and get your data faster.

SELECT Head 8/50 µL

Order Number: OL3318-26-410

Plate format: 96, 384 wells and tubes

Volume range: 1 – 50 µL

Specified volume range precision:

3 – 5 µL ≤ 3%

5 – 50 µL ≤ 2%

SELECT Head 8/250 µL

Order Number: OL3318-26-420

Plate format: 96, 384 wells and tubes

Volume range: 5 – 250 µL

Specified volume range precision:

10 – 25 µL ≤ 3%

25 – 250 µL ≤ 2%

SELECT Head 8/1000 µL

Order Number: OL3318-26-430

Plate format: 96, 384 wells and tubes

Volume range: 10 – 1000 µL

Specified volume range precision:

25 – 100 µL ≤ 3%

100 – 1000 µL ≤ 2%



Explore the CyBio FeliX SELECT Head



© 2024 Analytik Jena GmbH+Co. KG
889-80001-2-F
Version 1.1 - en - 04/2025

Pictures: Analytik Jena GmbH+Co. KG

Subjects to changes in design and scope of delivery as well as further technical development.

Analytik Jena GmbH+Co. KG
Konrad-Zuse-Strasse 1
07745 Jena · Germany

Tel +49 3641 77 70
Fax +49 3641 77 9279
info@analytik-jena.com
www.analytik-jena.com

analytikjena
An Endress+Hauser Company