

PULSEspencer R

Technical Data



Technical Data

PULSEspencer R

General

The **PULSEspencer R** is engineered to deliver ultra-precise, contactless dispensing of reagents in volumes ranging from picoliters to microliters. Designed for life science researchers, it simplifies complex liquid handling tasks such as dose-response assays, qPCR setups, and compound screening. With its unmatched accuracy and speed, the **PULSEspencer R** supports miniaturized workflows, reduces reagent consumption, and accelerates assay development – all within a compact and intuitive platform.

PULSEspencer R (875-300010)

Weight	16.1 kg (35.4 lb)
Dimensions (L x W x H)	470 mm x 380 mm x 230 mm
Recommended footprint (incl. laptop)	850 mm x 380 mm x 230 mm
Dispensing Technology	<u>Non-contact</u> by applying thermally generated pressure pulses
Volume Range	8 pL – 200 µL* (approx. 7.4 orders of magnitude)
Precision	< 8%**
Dispensing Speed	10 nL across 96-well format in 20 sec
Compatible labware	12- to 1536-well plates (6-47 mm height)***, glass slides and custom labware
Supported Fluids	DMSO, aqueous solutions, proteomics solutions, PCR master mixes

* depending on type of fluid and dispensehead cassette

** for volumes >100 pL

*** plate limits due to evaporation: DMSO dispensing max. 32 plates, aqueous dispensing max. 8 plates or 1,536 wells
(means : 1x 1536-well plate or 4x 384-well plate or 8x 96-well plate)

Technical Data

PULSEspencer R

Physical data

Interface	■ USB port ■ RS-232 not documented
Fuses	2 × 5 × 20 mm T 4 A 250 V AC
Line voltage	100 – 240 V AC
Line Frequency	50/ 60 Hz
Power consumption	Max. 120 watts
Noise emission	< 75 db (A)
Operation conditions	■ Indoor use only ■ Permissible ambient temperature +20 °C to +30 °C ■ Permissible relative humidity 30% to 80%, non-condensing ■ Maximum altitude 3000 m ■ Pollution degree 2

Analytik Jena Dispenser Software

Specifications	Description
Operating system	■ Windows 11 ■ Windows 10 ■ Windows 8
	Minimum screen resolution of 1280 x 720 pixels
Memory requirements	256 GB drive (minimum) recommended
Memory (RAM)	8 GB RAM recommended
Interfaces	USB
Software language	English

Consumables

Dispensing cassettes for PULSEspencer R				
Product Name	RR8	RR8s	RR4	RR4s
Usage type	Single-use	Single-use	Single-use	Single-use
Description	8 dispenseheads of 20 µL total volume each		4 dispenseheads of 200 µL total volume each	
Qualified fluid classes	Aqueous solutions	DMSO, master mix, aqueous solutions <u>requiring</u> surfactant	Aqueous solutions	DMSO, master mix, aqueous solutions <u>requiring</u> surfactant
Order number	875-300725	875-300720	875-300745	875-300740
Packaging unit(s)	20	20	20	20
Warranty	1 year	1 year	1 year	1 year
Dead volume	0.5 µL	0.5 µL	0.5 µL	0.5 µL
Dispense volume/ Fill volume	all wells > 1 nL 4 – 20 µL all wells < 1 nL 4 – 10 µL	all wells > 1 nL 2 – 20 µL all wells < 1 nL 2 – 10 µL	all wells > 1 nL 50 – 200 µL	all wells > 1 nL 50 – 200 µL
CV	< 8%	< 8%	< 8%	< 8%
Dispensing volume range	8 pL – 20 µL	11 pL – 20 µL	1 nL – 200 µL	1 nL – 200 µL

This document is true and correct at the time of publication; the information within is subject to change. Other documents may supersede this document, including technical modifications and corrections.

© Analytik Jena GmbH+Co. KG